Altium

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

How Does Bleed Impact GC/MS Data and How Can It Be Controlled?

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



Brewing Excellence: Quantitating Over 200 Pesticides in Black Tea with Steady Performance and Maximized Uptime by GC/MS/MS

Authors
Anastasia A. Andrianova and
Limian Zhao
Agilent Technologies, Inc.

Abstract

This application note presents results for the sensitive and robust quantitation of 245 pesticides in black tea extract with the Agilent 7010D Triple Quadrupole Mass Spectrometer (GC/TQ) featuring a second-generation High Efficiency Source 2.0 (HES 2.0), that addresses the challenges posed by residual pesticide analysis in complex matrices. By optimizing sample preparation and using state-of-the-art GC/MS hardware, including ion source technology and midcolumn backflushing, excellent calibration performance and sensitivity at low-ppb levels were achieved. The method demonstrated exceptional ruggedness and robustness over 800 consecutive injections of a black tea extract spiked with pesticides at 2 ppb, with high precision and low RSDs, ensuring prolonged instrument uptime and maximum throughput. The demonstrated limits of quantitation (LOQs) were as low as 0.01 ppb for over a third of evaluated compounds, and the calibration range spanned up to five orders of magnitude while meeting SANTE 11312/2021 guidelines. This application note highlights intelligent GC/TQ features, such as early maintenance feedback and an instrument health status dashboard, maintaining confidence in the results for a high-throughput analysis. The updated data acquisition platform provides enhanced user experience, including a new implementation of the retention time locking functionality.

Author
Brooke C. Reaser
Agilent Technologies, Inc.

Enhanced Longevity and Revolutionized Robustness for the Sensitive Detection of 190 Pesticides over 800 Injections

Abstract

Multiresidue pesticide analysis has become one of the most difficult but important analytical challenges for those using gas chromatography and mass spectrometry. The Agilent 8890 gas chromatograph (GC) coupled with the Agilent 7010 triple quadrupole mass spectrometer (GC/TQ) with a high efficiency source 2.0 (HES 2.0) upgrade is an analytically accurate, robust, and reproducible instrument for multiresidue pesticide analysis of complex samples. The analysis of 190 pesticides in a spinach extract was conducted using an Agilent QuEChERS extraction kit across 800 injections. Only GC inlet maintenance was required over the duration of the injections. The instrument configuration that enabled robust performance included a multimode inlet, a mid-column backflushing configuration, and the HES 2.0. No degradation of the analytical method, sensitivity, or instrument performance occurred, allowing for the high-throughput, accurate, robust, and sensitive detection of pesticides in spinach.

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