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LC - Acquisition

File Home Control

Take Release Status Method Single Sample Sequence Ad Hoc Mode Copy Delete Run Queue Activity Log Instrument Status Online Signals Solvent Monitor Spectrum Sequence Creation Template Status Windows

Sequence - Sequence1.sqx

General Properties Run Options Injections Table

#	Action	Sample type
1	☒ Inject	Blank
2	☒ Inject	Cal. Std.
3	☒ Inject	Cal. Std.
4	☒ Inject	Sample
5	☒ Inject	Sample
6	☒ Inject	Sample
7	☒ Inject	Sample
8	☒ Inject	Sample

Sample custom parameters Compound

* Indicates Mandatory Entry

Parameter Value

Compartment 1

1.00 L
Actual Volume

Purging bottle

Consumption (est.): 0.00 L
Leftover volume (est.): 1.00 L

No Bottle
Actual Volume

Consumption (est.): 0.00 L
Leftover volume (est.): 1.00 L

No Bottle
Actual Volume

Consumption (est.): 0.00 L
Leftover volume (est.): 1.00 L

No Bottle
Actual Volume

Consumption (est.): 0.00 L
Leftover volume (est.): 1.00 L

0.15 L
Actual Volume

Acetonitrile

Consumption (est.): 0.00 L
Leftover volume (est.): 0.15 L

1.00 L
Actual Volume


Water

Consumption (est.): 0.00 L
Leftover volume (est.): 1.00 L


Result path: C:\CDSP\Projects\SLS\Results
Result name:
Save result as: One result set

Validate Run

Current user: SYSTEM (SYSTEM), Active project: SLS

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
InfinityLab Assist Troubleshooting

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Not Ready

Health

**Maintenance**
Maintenance guides and tools for common tasks.

**Diagnostics**
Diagnostic tasks and wizards.

**Insights**
EMF counters and statistics about the instrument.

**Troubleshooting**
Guides and information for when the instrument runs into a problem.

**Log**
Instrument actions are recorded here.

Connection to Assist Hub is lost, trying to reconnect...

← Replace Deuterium Lamp

Overview

Results

VWD
G7114A:PH11287092

Time Required80 - 90 min

Steps2

Tools & Parts Required

Safety

WARNINGCAUTION

This procedure describes the replacement of the deuterium lamp.

Start Maintenance

Error

←

Assisted Troubleshooting: Select Cause

Most probable cause

Clogged needle or needle seat

Other causes

Clogged column

Misconfigured method parameters

Clogged heat exchanger

Clogged capillary

Cancel

Notifications

?

×

Binary Pump G7132A:LPBBV00014

less than a minute ago

^

▲

✖

Pressure above high pressure limit (600 bar)

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
InfinityLab Certified ID Vials

Sample Linking Station
with Sample Linking Software

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Analyzing Iohexol by Compendial Method Produces Excellent Precision with Shallow Gradients Using an Agilent 1260 Infinity III LC

Abstract
The separation of Iohexol and its impurities is a process performed for the analysis of pharmaceuticals. Several of these methods require conditions such as shallow gradients to achieve the required resolution and precision. This paper describes a method for the separation of Iohexol and its impurities using a shallow gradient method. The method involves the use of a compendial method to achieve the required resolution and precision. The method involves the use of a compendial method to achieve the required resolution and precision. The method involves the use of a compendial method to achieve the required resolution and precision.

Gradient	0 min	1% B
	60 min	13% B
	60.01 min	99% B


Figure 2. Chromatogram overlay (N = 8) of the separation using the compendial method.

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