UV-Vis-(NIR) spektrometrie

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
CARY 60 APLIKACE

Application Note
BioPharma/Pharma

Concentration Analysis of Ultra-microvolume Samples by UV-Vis Spectroscopy

Simple, precise, and nondestructive concentration analysis with Agilent Cary 60 UV-Vis and TrayCell 2.0

Authors

Geethika Weragoda,
Wesam Alwan, and
Fabian Zieschang
Agilent Technologies, Inc.

Abstract

The Agilent Cary 60 UV-Vis spectrophotometer fitted with a TrayCell 2.0 Ultra-Microvolume Cell was used for the analysis of protein and nucleic acid samples. The non-destructive measurement of microliter-sized sample volumes is beneficial in terms of sample preservation, sample handling, and eliminating the need to perform tedious and error-prone dilutions. Using the TrayCell 2.0 is quick and easy, resulting in workflow improvements for users.

Application Note
Environmental

SPADNS Colorimetric Method to Quantify Fluorides in Water

Convenient and precise concentration measurements with the Agilent Cary 60 UV-Vis spectrophotometer

Authors

Rogelio García and
Geethika Weragoda
Agilent Technologies, Inc.

Abstract

The Agilent Cary 60 UV-Vis spectrophotometer was used to quantify fluoride in water as described in the SPADNS colorimetric method. The Concentration module in Agilent Cary WinUV software was used for single wavelength absorbance measurements and data analysis, reducing time-consuming data workup procedures. The calibration curve has a linear analytical range between 0 to 1.4 mg/L, which can be used in analyzing water samples.

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





Benefits of the Cary 3500 Multicell UV-Vis for Protein Analysis

Productivity and reproducibility improvements for qualitative and quantitative measurements of very low volumes

Authors
Kevin Grant and Matt Quinn
Agilent Technologies,
Australia

Introduction

Measurements made by a UV-Vis spectrophotometer provide a fast quality control check of a sample. The measurements can also be used to estimate specific activity or estimate yield after purification, and to identify fractions containing protein or amino acids.

Proteins containing amino acids with aromatic side chains absorb light at approximately 280 nm. This fact enables quantitative analysis by a UV-Vis spectrophotometer. Following Beer-Lambert's law ($I = \epsilon \cdot c \cdot l$), the concentration of these amino acids can be determined when the absorbance is known. Often a protein sample is limited, and being able to do these measurements at a minimal volume is important for sample preservation. Performing a scan can provide information about potential contaminants.

The Cary 3500 Multicell UV-Vis has an integrated, permanently aligned sample holder, ideal for reproducible and reliable qualitative, and quantitative measurements with ultra-microvolume cuvettes (Figure 1).

Results

The hydrolysis of pNPA involves the removal of the acetate group under basic conditions. When the conditions are set such that water is in excess, then the reaction can be considered as pseudo first order. As the pH of the PBS buffer was set to 7, a slower reaction rate can be expected, and the second order behavior dominates.

Temperature effect

Wavelength scans of pNPA hydrolysis at four different temperatures are shown in Figure 2. There is a clear relationship between the temperature of the sample and the reaction rate as measured by PNP production (Figure 2). At 80 °C the isosbestic points are evident, indicating a direct conversion from pNPA to PNP. The wavelength scanning kinetics data for the 80 °C experiment was used to produce an absorbance versus time plot, using the PNP peak at 408 nm (Figure 3). The second order rate calculation (built in to Cary UV Workstation software) was then used to determine the second order rate constant (k) for the reaction as $k = 883.194 \text{ (1/min.mol)}$.

Wavelength scans

By performing wavelength scans over time, the consumption of pNPA and PNP production can be observed (Figure 2). The entire wavelength range provides additional information that might otherwise be missed if only a single wavelength was analyzed. For example, the possible presence of reaction intermediates, as well as subtle changes in the sample, and detection of isosbestic points as shown in Figure 2.


Figure 2. Wavelength scans over time over the wavelength range 220 - 520 nm, collected for 30 minutes after the reaction was initiated by mixing the two reagents. Top left is at 20 °C, top right is 40 °C, bottom left is 60 °C and bottom right is 80 °C.

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