

Characterization of Monoclonal Antibody (mAb) using Shimadzu Q-TOF LCMS 9030

Shannie Tay ¹, Max Kosok¹, Yu Jie Lee²
(1) Shimadzu (Asia Pacific), Singapore (2) National University of Singapore, Singapore.

1. Introduction

- ◆ Monoclonal antibody (mAb) (Fig.1) drugs are gaining traction in the biopharmaceutical industry. However, due to their complexity, they often require thorough characterization to assess Critical Quality Attributes (CQAs) during development stages and for regulatory approval. Typically, High Resolution Mass Spectrometry (HRMS) Quadrupole Time-of-Flight (Q-TOF) LCMS-9030 (Fig.2) is employed for mAbs characterization.
- ◆ During process development, mAbs undergo modifications like oxidation and deamination, impacting their safety and efficacy. Hence, post-translation modification (PTM) studies are crucial. Peptide mapping, a CQA method, confirms the amino acid sequence and identifies modification sites, aiding in assessing safety and efficacy..
- ◆ Intact and subunit analysis are conducted to ascertain the molecular weight and overall structure of the mAb. Additionally, peptide mapping is employed to achieve a thorough understanding of the sequence coverage of the mAb, while simultaneously identifying and locating post-translational modification sites.




Fig. 1 Structure of mAb



Fig. 2 QTOF-9030

2. Methods

In this study, the designated reference mAb used was the NISTmAb reference material 8671, a humanized IgG1k monoclonal antibody. For subunit analysis, the Rapid™ PNGase F enzyme was employed to effectively remove glycans from the heavy chain. For sequence confirmation and PTM analysis, the mAb underwent enzymatic digestion using trypsin to generate smaller peptide fragments. Analysis of the samples was performed using Shimadzu Q-TOF LCMS-9030 instrument, and Protein Metrics was employed for subsequent data processing.

3. Results




Fig. 3 MS Spectra of modified peptide with loss of amino acid Lysine (K) at heavy chain amino acid position 450




Fig. 4 MS Spectra of unmodified wildtype peptide without loss of amino acid Lysine (K) at heavy chain amino acid position 450

Protein name	Start AA	Var. Pos. Protein	Sequence (unformatted) †	Mod. Names	Mod. AAs	MS Id	1
						MS Alias name	NIST (%)
Heavy Chain	1	1	QVTLR	Gln->pyro-Glu/-17.0265	Q		99
			QVTLRESGPAIVKPTQTTLTCTFSGFSLTAGMSVGWIR	Gln->pyro-Glu/-17.0265	Q		1.04
	162	162	STSGGTAALGCLVKDYFPEPTVSWNSGALTSQVHTFPAVLQSSGLYSLSS	Deamidated/0.9840	N		2.81
	162,178	162,178	STSGGTAALGCLVKDYFPEPTVSWNSGALTSQVHTFPAVLQSSGLYSLSS	Deamidated/0.9840	N,Q		12.4
	252	255	DTLMISR	Oxidation/5.9949	M		2.88
	278	289	FNWYVDGVEVHNAK	Deamidated/0.9840	N		2.5
			TKPREEQYNSTYR	NGlycan/1079.4017	N		5.57
				NGlycan/1241.4545	N		23.2
				NGlycan/1403.5073	N		9.23
				NGlycan/1444.5339	N		33.4
				NGlycan/1506.5867	N		17.2
				NGlycan/1768.6395	N		3.86
				NGlycan/876.3223	N		1.9
				NGlycan/1038.3751	N		5.61
				NGlycan/1241.4545	N		189
	296	300	EEQYNSTYR	Deamidated/0.9840	Q		0.659
		314	VVSVLTVLHQDWLNGK	Deamidated/0.9840	Q		2.1
			VVSVLTVLHQDWLNGKEYCK	Deamidated/0.9840	N		0.942
	305	318	VVSVLTVLHQDWLNGKEY	Deamidated/0.9840	N		26
			VVSVLTVLHQDWLNGKEYCK	Deamidated/0.9840	N		2.2
	326	328	VSNKALPAIEK	Deamidated/0.9840	N		0.767
		348	EPQVYTLPPSREEMTK	Deamidated/0.9840	E		1.79
	348		EPQVYTLPPSREEMTKNQVSLTCLVK	Glu->pyro-Glu/-18.0106	E		1.59
		364	EPQVYTLPPSREEMTKNQVSLTCLVK	Deamidated/0.9840	N		2.1
		364	NQVSLTCLVK	Deamidated/0.9840	N		3.14
		364	NQVSLTCLVK	Deamidated/0.9840	Q		3.14
Light Chain			GFYPSDIAVEWESNGQPENNYK	Deamidated/0.9840	N		13.8
		387	GFYPSDIAVEWESNGQPENNYK	Deamidated/0.9840	N		4.95
			GFYPSDIAVEWESNGQPENNYK	Deamidated/0.9840	N,Q		6.51
		387,389	GFYPSDIAVEWESNGQPENNYKTPPVLDSGGSFFLYSK	Deamidated/0.9840	N,Q		5.57
			GFYPSDIAVEWESNGQPENNYK	Deamidated/0.9840	Q		11.7
		389	GFYPSDIAVEWESNGQPENNYKTPPVLDSGGSFFLYSK	Deamidated/0.9840	Q		3.8
		389,392	GFYPSDIAVEWESNGQPENNYK	Deamidated/0.9840	Q,N		6.51
			GFYPSDIAVEWESNGQPENNYK	Deamidated/0.9840	N		11.7
		392	GFYPSDIAVEWESNGQPENNYKTPPVLDSGGSFFLYSK	Deamidated/0.9840	N		3.8
	413	437	LTVDKSRWQGNVFCSCVMHEALHHYTK	Deamidated/0.9840	N		3.1
	420	437	WQQGNVFCSCVMHEALHHYTK	Deamidated/0.9840	N		5.4
	443	450	SLSLSPGK	Lys-loss/-128.0950	K		78.4
	108	136	TVAAPSVFIFPPSDEQLKSGTASVCLNNFYPR	Deamidated/0.9840	N		4.26
	126	136	SGTASVCLNNFYPR	Deamidated/0.9840	N		5.16
		151,154	VQWKVDNALQSGNSQESVTEQDSKDSTYSLSSTLTLSK	Deamidated/0.9840	N,Q		21.3
		157	VQWKVDNALQSGNSQESVTEQDSKDSTYSLSSTLTLSK	Deamidated/0.9840	N		4.03
	207	209	SFNRGEC	Deamidated/0.9840	N		2.33

Table 1 Amino acid site relative abundance %

- ◆ **Table 1** gives the summary of the overall calculated amino acid site relative abundance percentage.
- ◆ Protein Metrics software automatically calculates the relative percentage of PTM of each individual modified amino acid on the peptides based on the ratio of extracted ion chromatograms (XICs) corresponding to the modified and unmodified wildtype peptides.




Fig.5 Deconvoluted spectra of intact NISTmAb with modifiers




Fig.6 & 7 Deconvoluted spectra of NISTmAb light and heavy chain with modifiers

		Sample name	
		Peak #	1
Name	Mass	Expected mass	Delta mass from calculated
Antibody , G2F/G2F	148699.234	148685.875	13.36
Antibody , G1F/G2F + 4 Oxi + Glycation	148745.795	148749.775	-3.98
Antibody , G1F/G2F	148527.284	148523.775	3.51
Antibody , G1F/G1F + 6 Oxi	148471.338	148457.575	13.76
Antibody , G1F/G1F + 4 Oxi + Glycation	148582.439	148587.575	-5.14
Antibody , G1F/G1F + 3 Oxi	148416.051	148409.575	6.48
Antibody , G1F/G1F	148355.508	148361.575	-6.07

Table 2 Summary of NISTmAb intact mass analysis of theoretical vs. observed mass

		Sample name	NIST Subunit Analysis	NIST Subunit Analysis
		Peak #	1	2
Name	Mass	Expected mass	Delta mass from calculated	
Light Chain , Reference	23128.974	23123.593	5.38	
Light Chain , pyroQn	23097.992	23105.593	-7.6	
Light Chain , + 3 Oxi	23167.58	23171.593	-4.01	
Heavy Chain , Hexose	49613.982	49614.645		-0.66
Heavy Chain , + 6 Oxi	49544.999	49548.645		-3.65
Heavy Chain , + 3 Oxi	49501.9	49500.645		1.25

Table 3 Summary of NISTmAb subunit analysis of theoretical vs. observed mass

- ◆ Protein Metrics software automatically deconvolutes the m/z to give accurate molecular weight measurement of the intact mass and subunit analysis of the mAb, based on the default disulfide bond structure. By using high-resolution mass spectrometry, clear and distinguishable peak spectra can be obtained, which allows the protein metrics software to easily calculate the accurate molecular weight of the proteins with the glycoforms and modifiers (Fig. 5 to 7).

- ◆ **Tables 2 and 3** provide a comprehensive summary of the NISTmAb measured mass compared to the theoretical delta mass differences, which is generated by protein metrics software. The data demonstrates the overall calculated delta mass difference is less than 20 delta ensuring accurate characterization of the mAb.

4. Conclusion

- ◆ During drug development, it's vital to monitor alterations in PTM amino acid site relative abundance percentages, like methionine (M) oxidation or asparagine (N) deamidation, as they can affect mAbs' efficacy, stability, and immunogenicity.

- ◆ Intact mass and subunit analysis determine the average mass of mAb, offering insights into their structural integrity crucial for drug product development. Subunit analysis assesses the light and heavy chain molecular weight following intra-chain disulfide bond reduction.

- ◆ The Shimadzu Q-TOF LCMS-9030 instrument is highly effective in generating precise data, which is crucial for the thorough characterization of mAb.

Reference

- 1) Monoclonal Antibody (mAb) Intact Mass and Subunit Analysis Using Shimadzu Q-TOF LCMS-9030, 04-AD-0296-EN
- 2) A Fast and Simple Workflow for Monoclonal Antibody (mAb) Post-Translation Modifications (PTM) Study Using Shimadzu LCMS-9030 Q-TOF, AD-0295-EN

Disclaimer:

The products and applications in this presentation are intended for Research Use Only (RUO). Not for use in diagnostic procedures.

***** All authors declare that they have no conflicts of interest. *****