

Simplifying EPA 1633 analyses with an improved dual bed solid-phase extraction method

Jason Hoisington¹, Alexis Shelow¹, Jason Thomas¹, Colton Myers¹, Diego A. Lopez¹

¹Restek Corporation, Bellefonte, PA

PFAS and EPA 1633


Per- and Polyfluoroalkyl Substances (PFAS) are synthetic organofluorine chemicals that have been linked to a range of health effects, including decrease immune response, vaccine effectiveness, cancer risk, among others.

EPA Method 1633 was finalized early last year, and it represents a comprehensive analysis of 40 target PFAS compounds in different aqueous, solids, and tissue matrixes. This method uses a weak anion exchange (WAX) solid phase extraction (SPE), coupled with a dispersive graphitized carbon black (GCB) clean-up step prior to LC-MS/MS analysis on a C18 column.

Statement of Problem


1. Dispersive GCB can be a cumbersome and laborious step that requires further filtration before LC-M/MS analysis.
2. EPA Method 1633 requires analysis of large volume (up to 500 mL) samples of complex aqueous matrices with a varying degree of suspended solids that tend to clog the SPE cartridge, even with addition of glass wool or settling and centrifugation steps. A secondary SPE tube is necessary if clogging occurs.

Figure 1. Inter-cartridge variation from the amount of wool used to fill up to “half the height” of the SPE tube. Manual packing carried out in triplicates by four different lab technicians.



Solution Approach


A dual bed WAX/GCB cartridge that leverages an integrated Filter Aid; therefore, eliminating the need for dispersive carbon clean up, as well as the requirement to hand-pack deactivated glass wool in each tube.



Sample Preparation Workflow

Samples for precision, accuracy, and method detection limit (MDL) studies were prepared in polypropylene bottles using 500 mL deionized water spiked with 25 μ L of extracted internal standards (EIS, Wellington Laboratories, Cat. # MPFAC-HIF-ES) as per EPA Method 1633, Section 11.2.4. Four samples for precision and recovery analysis were spiked with 200 μ L of native PFAS standards (Wellington Laboratories, Cat. # EPA-1633STK), giving pre-extraction concentrations of 100-2500 ng/L. Seven MDL samples were spiked with 20 μ L of a 20:1 dilution of the native standard, giving pre-extraction concentrations of 0.5-12.5 ng/L. Seven blank replicates were prepared for the MDL study as well. The MDL samples were prepared and analyzed over three days.


Figure 2. Sample preparation procedure per EPA 1633 using Thermo Scientific AutoTrace 280 PFAS System.



LC-MS/MS Parameters

Figure 3. Analytical conditions (Waters Xevo TQ-S with ACQUITY Premier UPLC) and chromatogram showing all target PFAS and isotopically labeled compounds.

Analytical Column	1.8 µm Force C18 50x2.1 mm id	
PFAS Delay Column	5 µm PFAS Delay 50x2.1 mm id	
Mobile Phase A	Water, 5mM Ammonium Acetate	
Mobile Phase B	Methanol	
Gradient	Time (min)	%B
	0.00	20
	6.00	95
	6.60	95
	6.61	20
	7.50	20
Flow Rate	0.4 mL/min	
Injection Volume	3 µL	
Column Temp.	40°C	
Ion Mode	Scheduled MRM in negative ESI	



Results and Discussion

Table 1. Results from MDL, Blank, Accuracy, and Precision experiments for native PFAS.

Compound	Abbreviation	MDL (g/mol)	Blank (%)	Accuracy (%)	%RSD
Perfluorooctanoic acid	PFBA	0.34	ND	112	5
Perfluorooctanoic acid	PFBA	0.34	ND	112	5
3-Perfluorooctanoyl propionic acid	PFMPA	0.20	ND	109	6
3-Perfluorooctanoyl propionic acid	3-PFCTA	0.31	ND	95	6
Perfluorooctanoic acid	PFBA	0.28	ND	94	6
Perfluorobutanoic acid	PFBS	0.19	ND	101	5
Perfluoro-4-methylbutanoic acid	PFMBa	0.16	ND	111	6
Perfluoro-2-ethylhexanoic acid	PFHEa	0.15	ND	94	6
Nonfluoro-3,6-dioxahexanoic acid	NFDHa	0.48	ND	121	5
11H,12H,2H-perfluorooctanoic acid	4-2 FTS	0.36	ND	104	5
2H,3H,4H,3H-perfluorooctanoic acid	5-5 FCTA	1.48	ND	102	6
Perfluorooctanoic acid	PFBA	0.38	ND	113	5
Perfluorooctanoic acid	PFBS	0.07	ND	122	4
Hexafluorocyclopentene oxide dimer acid	HFPO-DA	0.40	ND	114	12
Perfluorooctanoic acid	PFBA	0.36	ND	109	6
Perfluorooctanoic acid	PFHS	0.07	ND	104	7
4,8-Dioxia-3H-perfluorooctanoic acid	ADONA	0.39	ND	99	9
11H,12H,2H-perfluorooctanoic acid	4-2 FTS	0.36	ND	107	4
3-Perfluorooctyl propionic acid	3-PFCTA	2.32	ND	100	6
Perfluorooctanoic acid	PFHS	0.76	ND	121	8
Perfluorooctanoic acid	PFDA	0.32	ND	104	2
Perfluorooctanoic acid	PFDA	0.20	ND	104	2
Perfluorooctanoic acid	PFOA	0.28	ND	109	6
Perfluorooctanoic acid	PFOA	0.28	ND	109	6
Perfluorooctanoic acid	PFOA	0.28	ND	109	6
Perfluorooctanoic acid	PFOA	0.28	ND	109	6
Perfluorooctanoic acid	PFOA	0.28	ND	109	6
N-methyl perfluorooctanoic acid	NMFOA	0.67	ND	104	5
9-Chlorododecafluoro-3-oxanone-1-sulfonic acid	9-Cl-PFDSO	0.61	ND	85	11
Perfluorooctanoic acid	PFBS	0.25	ND	79	7
Perfluorooctanoic acid	PFBS	0.18	ND	117	8
N-methyl perfluorooctanoic acid	NMFOA	0.14	ND	112	11
11H,12H,2H-perfluorooctanoic acid	4-2 FTS	0.36	ND	113	11
Perfluorooctanoic acid	PFBA	0.26	ND	125	12
N-methyl perfluorooctanoic acid	NMFOA	0.15	ND	95	12
N-methyl perfluorooctanoic acid	NMFOA	0.23	ND	93	12
N-methyl perfluorooctanoic acid	NMFOA	0.23	ND	93	12
Perfluorooctanoic acid	PFDS	0.25	ND	94	23
11-Chlorododecafluoro-3-oxadecane-1-sulfonic acid	11-Cl-PDS	0.56	ND	77	25
Perfluorooctanoic acid	PFDA	0.27	ND	118	15
N-methyl perfluorooctanoic acid	NMFOA	1.11	ND	106	9
N-methyl perfluorooctanoic acid	NMFOA	1.10	ND	111	9
Perfluorooctanoic acid	PFDA	0.33	ND	108	15
Perfluorooctanoic acid	PFDS	0.13	ND	109	12
Perfluorooctanoic acid	PFDS	0.27	ND	109	12

Table 2. Results from Precision and Accuracy experiments for isotope dilution standards.

Compound	Abbreviation	Accuracy (%)	% RSD
Perfluoro-1-(1,2,3,4-13C ₄)butanoic acid	13C ₄ -PFBA	88	4
Perfluoro-1-(1,2,3,4,5-13C ₅)pentanoic acid	13C ₅ -PFPeA	66	5
Perfluoro-1-(1,2,3,4,5,6-13C ₆)hexanoic acid	13C ₆ -PFHx	97	3
Sodium 1H,1H,2H,2H-perfluoro-1-(1,2-13C ₂)hexane sulfonate	13C ₂ -d ₂ -FTS	95	8
Perfluoro-1-(1,2,3,4,6-13C ₅)hexanoic acid	13C ₅ -PFHxA	122	2
2,3,3,3-Tetrafluoro-1-(1,2,2,3,3,3-heptafluoropropyl)3-oxopropanoic acid	13C ₃ -PFPOD-AA	92	4
Perfluoro-1-(1,2,3,4-13C ₄)heptanoic acid	13C ₄ -PFHpA	107	4
Sodium perfluoro-1-(1,2,3-13C ₃)hexafluoroisobutanoate	13C ₃ -PFHIS	102	4
Sodium 1H,1H,2H,2H-perfluoro-1-(1,2-13C ₂)octanoate sulfonate	13C ₂ -d ₂ -FTS	83	5
Perfluoro-1-(13C ₃)nonanoic acid	13C ₃ -PFDA	106	5
Sodium perfluoro-1-(13C ₃)decanoate sulfonate	13C ₃ -PFDS	77	7
Perfluoro-1-(13C ₃)tridecanoic sulfonamide	13C ₃ -PFMTSA	94	11
N-methyl-D-3-perfluoro-1-octanesulfonamide;dicarboxic acid	D ₃ -NMFOSA	85	9
N-methyl-D-3-perfluoro-1-octanesulfonamide;dicarboxic acid	D ₅ -NEFOSA	87	11
Perfluoro-1-(13C ₃)nonanoic acid	13C ₃ -PFDA	83	7
Perfluoro-1-(1,2,3,4,5,6-13C ₆)undecanoic acid	13C ₆ -PFUA	80	7
Sodium 1H,1H,2H,2H-perfluoro-1-(1,2-13C ₂)decanoate sulfonate	13C ₂ -d ₂ -FTS	104	5
N-methyl-D-7-perfluoro-octanesulfonamide;octanoic acid	D ₇ -NMFOSE	108	13
N-methyl-D-9-perfluoro-octanesulfonamide;octanoic acid	D ₉ -NEFOSE	88	17
N-methyl-D-3-perfluoro-1-octanesulfonamide	D ₃ -NMFOSA	84	8
N-methyl-D-5-perfluoro-1-octanesulfonamide	D ₅ -NEFOSA	111	10
Perfluoro-1-(1,2,3,4,5,6,7-13C ₇)tridecanoic acid	13C ₇ -PFUDA	78	13
Perfluoro-1-(1,2-13C ₂)dodecanoic acid	13C ₂ -PFDDA	45	24
Perfluoro-1-(1,2-13C ₂)tridecanoic acid	13C ₂ -PFDDA	36	21




Figure 4. Accuracy (%) for all target PFAS compounds from four replicate spikes. Accuracy ranged from 77% to 125% of the spiked value, meeting the requirements in Table 5 of EPA Method 1633.




Figure 5. Average volume of 500 mL wastewater extracted before SPE clogging. Comparison between Filter Aid, glass wool, and control. ASTM substitute wastewater (ASTM International, D5905-98[2018]) diluted to ~100 mg/mL suspended solids was used.

Summary

Nondrinking water matrices present unique challenges to labs using SPE sample preparation. In this study, a dual bed WAX/GCB cartridge format exhibited consistency with the acceptance criteria in EPA method 1633. Furthermore, the *Filter Aid* greatly mitigates the clogging risk by this challenging matrix, allowing for 100% extraction of the recommended sample volume under method guidelines.

COI Disclosure: Authors are employed by Restek Corporation.