

Sensitive and Rapid Analysis of Inorganic Arsenic in Food and Animal Feed Samples Using LC-ICP-MS

Kana Matsuno and Kosuke Naka
Shimadzu Corporation, Analytical & Measuring Instruments Division, Japan



1. Introduction

- ◆ Inorganic arsenic (iAs) is more toxic than organic arsenic, and its monitoring in food and feed is required worldwide.
- ◆ Because organic arsenic and chlorine in food matrices can interfere with iAs determination, effective chromatographic separation is essential. This study quantified iAs in various foods and feeds by LC-ICP-MS, following EN 16802:2016¹⁾ and EN 17374:2020²⁾.

2. Instrument Configuration

- ◆ Samples were analyzed using the LC-ICP-MS system combining the ICPMS-2040/2050 with a Nexera™ XS inert, controlled by LabSolutions™ ICPMS TRM software (Fig. 1). The column was HAMILTON PRP-X100 (10.0 μm, 250 mm × 4.1 mm), and the eluent was 100 mM ammonium carbonate + 3% v/v methanol + 400 mg/L sodium sulfate (pH 9.3).
- ◆ The mini torch for ICP-MS reduces argon gas consumption to 11 L/min, roughly two-thirds that of a conventional torch, lowering running costs in lengthy LC-ICP-MS analyses.



Fig. 1 Workflow of this study

3. Regulatory Values and Detection Limits

- ◆ The limit of detection (LOD) were sufficient even for the strictest iAs limit under the Baby Food Safety Act 2021 (Table 1).

Samples	Regulatory Values (mg/kg)				LODs in samples (mg/kg)	LOD in analytical solution (μg/L)
	(EU) 2023/915 ³⁾	(EU) 2025/1891 ⁴⁾	The Baby Food Safety Act of 2021 ⁵⁾	DIRECTIVE 2002/32/EC ⁶⁾		
Brown Rice	0.25	-	-	-	0.0004	0.02
Hijiki	-	-	-	-	0.04	
Rice Cereal	0.02	-	0.015	-	0.0004	
Apple Jelly	0.02	-	0.01	-	0.0004	
Strained Corn	0.02	-	0.01	-	0.0004	
Infant Formula	0.02	-	0.01	-	0.0004	
Sea Bream	-	0.1	-	-	0.0004	
Tuna	-	0.1	-	-	0.0004	
Squid	-	0.05	-	-	0.0004	
Ham	-	-	-	-	0.0004	
Orange Juice	0.02	-	-	-	0.0002 (mg/L)	
Japanese Sake	-	-	-	-	0.00004 (mg/L)	
Animal Feed	-	-	-	2	0.002	

4. Oxidation from As(III) to As(V)

- ◆ Without pretreatment, As(III) appeared as a broad peak, but after pretreatment a peak was obtained at the same retention time as As(V).
- ◆ When the 5 μg/L As(III) check standard (prepared from As(III)) was quantified using the iAs calibration curve prepared from As(V) standards, the measured concentration was 4.83 μg/L. These results confirm that As(III) is oxidized to As(V) by the pretreatment (Fig. 2).




Fig. 2 Chromatograms for As(III) Standard Solutions (shown with baseline offset)

5. Removal of Cl Interference

- ◆ The Cl and iAs peaks were resolved under the chromatographic conditions used.
- ◆ The use of helium collision mode effectively suppressed ⁴⁰Ar³⁵Cl interference at m/z 75 (Fig. 3).




Fig. 3 Removal of Interference derived from Cl

(a) Separation of iAs from Cl, (b) Removal of ⁴⁰Ar³⁵Cl Interference in He Collision

6. Separation of iAs from Organic As in Samples

- ◆ Although tuna and squid contained large amounts of organic arsenic (AsB), the analysis of iAs in the 5 μg/L spiked seafood samples was not affected (Fig. 4).




Fig. 4 Chromatograms of Seafood Spiked Sample
(a) Chromatograms of Seafood Spiked Sample, (b) Enlarged Chromatograms of Seafood Spiked Sample (shown with baseline offset)


7. Validation and Quantitative Analysis

- ◆ Quantified values for all three Certified Reference Materials (CRM) fell within their certified ranges (Table 2), confirming method validity.

Table 2 Results of CRM (iAs)			
Samples	LODs in samples (mg/kg)	Certified Values (mg/kg)	Quantitative Results (mg/kg)
Brown rice NMIJ 7532-a	0.0004	0.298 ± 0.008	0.291
Brown rice NMIJ 7533-a	0.0008	0.530 ± 0.016	0.527
Hijiki NMIJ 7405-b	0.04	24.4 ± 0.7* ¹	24.7

*1: Certified value as As(V)

- ◆ For the samples other than the CRMs, spike-recovery tests yielded 96–113% iAs recovery across food and feed matrices (Fig. 5), confirming method accuracy and reliability.
- ◆ Table 3 shows the quantified iAs results. Most samples contained little or no detectable iAs, while rice cereal and animal feed showed higher levels. This may be related to their grain-based ingredients.



8. Conclusion

- ◆ Quantitative analysis of inorganic arsenic using the Shimadzu LC-ICP-MS system was conducted on 14 different food and feed matrices.
- ◆ Inorganic arsenic was separated from organic arsenic and chlorine in various food and feed samples within 6.5 minutes.
- ◆ CRM analysis and spike recovery tests confirmed accurate inorganic arsenic analysis in a wide range of food and feed samples. Detection limits were low enough for even strict baby food regulations.

Reference

- 1) EN 16802:2016: Determination of inorganic arsenic in foodstuffs of marine and plant origin by anion-exchange HPLC-ICP-MS
- 2) EN 17374:2020: Determination of inorganic arsenic in animal feed by anion-exchange HPLC-ICP-MS
- 3) COMMISSION REGULATION (EU) 2023/915 of 25 April 2023 on maximum levels for certain contaminants in food and repealing Regulation (EC) No 1881/2006
- 4) Commission Regulation (EU) 2025/1891 of 17 September 2025 amending Regulation (EU) 2023/915 as regards maximum levels of inorganic arsenic in fish and other seafood
- 5) US House of Representatives, The Baby Food Safety Act of 2021, accessed Jan 2026
- 6) DIRECTIVE 2002/32/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 7 May 2002

Nexera and LabSolutions are trademarks of Shimadzu Corporation or its affiliated companies in Japan and/or other countries.